



***IP-HCMC Immunology course 2005***

***Practical course***

**Cytokine production determination by intracellular flow  
cytometry**

**February 1<sup>st</sup>, 2005**

Sylviane PIED, Institut Pasteur Paris

Daniel SCOTT-ALGARA, Institut Pasteur Paris

The purpose of this practical study is to determine the cytokine levels by intracellular flow cytometry method. This method can be used in both basic and clinical research for human and animal models.

The principle of this technique is to quantify cytokine-producing cells by detecting cytokines by antibody staining and flow cytometry analysis.

The term "flow cytometry" derives from the measurement (meter) of single cells (cyto) as they flow past a series of detectors. The fundamental concept is that cells flow one at a time through a region of interrogation where multiple biophysical properties of each cell can be measured at rates of over 1000 cells per second. These biophysical properties are then correlated with biological and biochemical properties of interest. The high through-put of cells allows for rare cells, which may have inherent or inducible differences, to be easily detected and identified from the remainder of the cell population.

Flow Cytometry involves the use of a beam of laser light projected through a liquid stream that contains cells, which when struck by the focussed light give out signals which are picked up by detectors. These signals are then converted for computer storage and data analysis, and can provide information about various cellular properties.

## Methodology

### I-Isolation of mouse splenocytes

- Mice are killed by exposition to CO<sub>2</sub> during at least 5 minutes
- Dissect the animal to take the spleen
- Dilacerate the spleen in sterile RPMI complemented with PS and transfer the cells in a 15 ml tube.
- After 5 minutes the cell suspension is decanted to a new 15 ml tube.
- Centrifuge 5 min at 300g and wash twice with RPMI
- Resuspend the pellet in 10 ml of complete RPMI culture medium (500 ml RPMI with glutamate supplemented with 50 ml fetal calf serum, 5 ml streptomycin/penicillin antibiotic solution, 5 ml Hepes 1M and 1 ml  $\beta$ 2 mercaptoethanol).
- Count cells on a hemocytometer.

### II-Cells stimulation

- **Note:** *The culture times indicated in this protocol are optimised for IL-10 and IFN $\gamma$ . Culture times for optimal staining of other cytokines will need to be determined empirically by the investigator depending on the type of cell used and the stimulation process.*
- Dilute cells to  $2.5 \times 10^6$  cells/ml into two separates culture flasks.
- Add 50 ng/ml PMA + 500  $\mu$ g/ml Ionomycin + 10  $\mu$ g/ml Brefeldin A to one flask.
- Incubate flasks at 37°C, 5%CO<sub>2</sub> for 3-4 hours.
- Recover cells after stimulation and wash 1X with PBS|| **WORK ON ICE!!!!**
- **Note:** *Brefeldin A is a fungal metabolite that disrupts the structure and function of the Golgi apparatus. Therefore, protein secretion is inhibited and newly synthesised proteins accumulate inside the cells. Ince it is toxic, Brefeldin A should not be added too long.*

- #### **IV-Intracytoplasmic cytokines staining**

## **Annexe: Reagents, solution disposable**

**Animals:** 8-10 weeks-old mice

### **Antibodies:**

- anti-CD4-PE (clone H129.19)
- anti-CD3-FITC (145-2C11)
- anti-IFN $\gamma$ -APC (clone XMG1.2)
- anti-IL-4-APC (clone 11B11)
- anti-Fc $\gamma$ II/III (clone 24G2)

All antibodies are from BD-Pharmingen.

### **Reagents:**

- Fetal calf serum
- Formaldehyde 4%
- PBS 1 X
- Ionomycin calcium salt (sigma ref. I-0634, 1mg)
- Brefeldin A (Sigma ref B7651, 5mg)
- Phorbol 12 Myristate 13 acetate (PMA, Sigma ref P-8139; 1mg)
- Saponin (Sigma ref S-7900)
- Paraphormaldehyde (PFA) 4%
- Washing and incubation solution: PBS 1X-3% FCS-0,1%NaN<sub>3</sub> (Azide).
- Fixation solution: PFA 2% final
- Permeabilisation buffer: Permash 1X (BD-Pharmingen)
- Or 1X PBS-3%FCS-0.5%Saponin

### **Disposables:**

- Absorbent paper
- Beakers
- Centrifuge tubes "falcon" 15 ml
- Facs Tubes
- 96 microwells plates
- Tissue culture flasks 75 cm<sup>3</sup>
- Gloves size L, M, S
- Micropipette tips 200-1000 $\mu$ l
- Micropipette tips 20-200 $\mu$ l
- Sterile pipettes 10 ml
- Sterile pipettes 5 ml
- Ice

### **Equipment:**

- Dissection tools
- Cell culture CO<sub>2</sub> incubator
- Cell culture Hood
- Flow cytometer

Malassez cytometer  
Micropipettes 1000µl  
Micropipettes 200 µl  
Micropipette 20 µl  
Microscope  
Multichannel pipette 200 µl  
Vortex  
Refrigerated centrifuge for 15 ml tubes